

Tomm7 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9POU1

Item Number: CM0007057



Introduction

Highly specific rabbit monoclonal TOMM7 antibody validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human and mouse samples. Ideal for mitochondrial import, PINK1/Parkin mitophagy, and TOM complex research.

[Learn More](#)

Parameter	Value
Product Name	TOMM7 Rabbit mAb
Remarks/Alias	TOM7; TOMM7
Species	Human
GeneID (Human)	54543
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 37-55 of human TOMM7 (NP_061932.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9POU1
Protein Weight	6kDa
Shipping	Ice bag