


Topoisomerase II Alpha/Top2A Rabbit Mab For Wb, If/Icc, If-P, Ihc-P, Elisa - P11388

Item Number: CM0011770



Introduction

Topoisomerase II alpha/TOP2A Rabbit monoclonal antibody for WB, IF/ICC, IHC-P, and ELISA. Validated in Human, Mouse, Rat. Immunogen: synthetic peptide from C-terminal region.

[Learn More](#)

Parameter	Value
Product Name	Topoisomerase II alpha/TOP2A Rabbit mAb
Remarks/Alias	TOP2; TP2A; TOP2A; TOP2alpha
Species	Human
Gene ID	7153
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1431-1531 of human Topoisomerase II alpha/TOP2A (NP_001058.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11388
Protein Weight	174kDa
Shipping	Ice bag