

Tp53Bp1 Rabbit Pab - Q12888

Item Number: CM0019150



Introduction

TP53BP1 Rabbit pAb validated for Western blot, IHC-P, and ELISA. Reacts with human, mouse, and rat. Recognizes the 214kDa DNA repair protein 53BP1, essential for double-strand break repair and class-switch recombination research.

[Learn More](#)

Parameter	Value
Product Name	TP53BP1 Rabbit pAb
Remarks/Alias	p202; 53BP1; TDRD30; p53BP1; TP53BP1
Species	Human
Gene ID (Human)	7158
Gene ID	7158
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 895-1031 of human TP53BP1 (NP_005648.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12888
Protein Weight	214kDa
Shipping	Ice bag