

Tp53I11 Rabbit Pab - O14683

Item Number: CM0017246

Introduction

Rabbit polyclonal antibody against human TP53I11 (PIG11) validated for WB, IHC-P, and ELISA, with cross-reactivity to mouse and rat. Designed for p53 signaling and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	TP53I11 Rabbit pAb
Remarks/Alias	PIG11; TP53I11
Species	Human
GeneID (Human)	9537
GeneID	9537
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-70 of human TP53I11 (NP_001245249.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14683
Protein Weight	21kDa
Shipping	Ice bag