

Tpm3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P06753

Item Number: CM0017038



Introduction

High-quality TPM3 Rabbit Polyclonal Antibody suitable for Western blot, immunofluorescence, and ELISA applications. Detects tropomyosin alpha-3 chain in human, mouse, and rat samples. Essential for cytoskeleton and muscle biology studies.

[Learn More](#)

Parameter	Value
Product Name	TPM3 Rabbit pAb
Remarks/Alias	TM3; TM5; TRK; CFTD; NEM1; TM-5; TM30; CAPM1; CMYP4A; CMYP4B; TM30nm; TPM3nu; TPMsk3; hscp30; HEL-189; HEL-S-82p; OK/SW-cl.5; TPM3
Species	Human
GeneID (Human)	7170
GeneID	7170
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-285 of human TPM3 (NP_689476.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06753
Protein Weight	33kDa
Shipping	Ice bag