


Traf4 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Buz4

Item Number: CM0029536



Introduction

High-quality rabbit anti-TRAF4 polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Detects endogenous TRAF4 in human, mouse, and rat samples. Ideal for studies on EGFR signaling, apoptosis, and cell differentiation.

[Learn More](#)

Parameter	Value
Product Name	TRAF4 Rabbit pAb
Remarks/Alias	CART1, MLN62, RNF83, TRAF4
Species	Human
GeneID (Human)	9618
GeneID	9618
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 251-470 of human TRAF4 (NP_004286.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BUZ4
Protein Weight	54kDa
Shipping	Ice bag