


Trail/Tnfsf10 Rabbit Mab - P50591

Item Number: CM0004811



Introduction

High-quality rabbit monoclonal antibody against human TRAIL/TNFSF10, validated for WB, IHC-P, IF/ICC, ELISA. Recognizes human, mouse, rat. Ideal for apoptosis and cancer research.

[Learn More](#)

Parameter	Specification
Product Name	TRAIL/TNFSF10 Rabbit mAb
Remarks/Alias	TL2; APO2L; CD253; TANCRC; TRAIL; Apo-2L; TNLG6A
Species	Human
GeneID (Human)	8743
GeneID	8743
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 114-281 of human TRAIL/TNFSF10 (NP_003801.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50591
Protein Weight	12 kDa/33 kDa
Shipping	Ice bag