


Trem11 Rabbit Polyclonal Antibody For Wb, Elisa - Q86Yw5

Item Number: CM0017971



Introduction

Rabbit polyclonal antibody against human TREML1 (TLT-1, triggered receptor expressed on myeloid cells-like protein 1) validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Suitable for innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	TREML1 Rabbit pAb
Remarks/Alias	TLT1; TLT-1; PRO3438; GLTL1825; dj238023.3; TREML1
Species	Human
Gene ID	340205
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 16-162 of human TREML1 (NP_835468.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q86YW5
Protein Weight	33kDa
Shipping	Ice bag