

Trim27 Rabbit Pab - P14373

Item Number: CM0024172

Introduction

TRIM27 rabbit polyclonal antibody validated for IHC-P and ELISA. Recognizes human, mouse, rat TRIM27, enabling research on ubiquitination, apoptosis, innate immunity, and autophagy pathways.

[Learn More](#)

Parameter	Value
Product Name	TRIM27 Rabbit pAb
Remarks/Alias	RFP; RNF76; TRIM27
Species	Human
Gene ID (Human)	5987
Gene ID	5987
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human TRIM27 (NP_006501.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14373
Protein Weight	58kDa
Shipping	Ice bag