

Trim47 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q96Ld4

Item Number: CM0019883

Introduction

Rabbit polyclonal antibody against human TRIM47, validated for WB, IF/ICC, ELISA; detects endogenous TRIM47 across human, mouse, rat. Suitable for ubiquitin pathway and cancer studies.

[Learn More](#)

Parameter	Specification
Product Name	TRIM47 Rabbit pAb
Remarks/Alias	GOA; RNF100; TRIM47
Species	Human
GenelD (human)	91107
GenelD	91107
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 356-530 of human TRIM47 (NP_258411.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96LD4
Protein Weight	70kDa
Shipping	Ice bag