


Trim5 Rabbit Pab For Wb, Elisa - Q9C035

Item Number: CM0020938



Introduction

High-quality TRIM5 Rabbit Polyclonal Antibody validated for WB and ELISA applications. Suitable for human and mouse samples. Recognizes endogenous TRIM5 (Q9C035) for retroviral restriction research.

[Learn More](#)

Parameter	Value
Product Name	TRIM5 Rabbit pAb
Remarks/Alias	RNF88; TRIM5alpha; TRIM5
Species	Human
Gene ID (Human)	85363
Gene ID	85363
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 100-200 of human TRIM5 (NP_149083.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9C035
Protein Weight	56kDa
Shipping	Ice bag