


Trim5 Rabbit Pab - Q9C035

Item Number: CM0021599



Introduction

TRIM5 Rabbit Polyclonal Antibody (pAb) validated for Western blot (WB) and ELISA. Detects endogenous TRIM5 in human, mouse, and rat samples; predicted molecular weight ~56 kDa. Suitable for antiviral innate immunity and autophagy research.

[Learn More](#)

Parameter	Value
Product Name	TRIM5 Rabbit pAb
Remarks/Alias	RNF88, TRIM5alpha, TRIM5
Species	Human
GeneID (Human)	85363
GeneID	85363
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 100-200 of human TRIM5 (NP_149083.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9C035
Protein Weight	56kDa
Shipping	Ice bag