

Trimethyl-Histone H3-K14 (Q16695 / P68431) Rabbit Polyclonal Antibody For Wb, If, Elisa

Item Number: CM0000343



Introduction

Rabbit polyclonal antibody specific for tri-methylated lysine 14 of Histone H3 (UniProt Q16695/P68431). Validated for WB, IF/ICC, ELISA, with reactivity in human, mouse, rat, and wide predicted cross-species coverage. Ideal for epigenetic and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	TriMethyl-Histone H3-K14 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; TriMethyl-Histone H3-K14
Species	Human
GeneID (Human)	H3
GeneID	8290/8350
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Histone H3 (NP_003520.1)
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag