

Trimethyl-Histone H3-K27 Rabbit Polyclonal Antibody (Q16695/P68431) For Wb, Chip, Elisa

Item Number: CM0000367



Introduction

Rabbit polyclonal antibody recognizing trimethylated lysine 27 of histone H3 (SWISS Q16695/P68431). Validated for WB, IP, IF/ICC, IHC-P, ChIP, ChIP-seq, ELISA. Suitable for human, mouse, rat, and other species.

[Learn More](#)

Parameter	Value
Product Name	TriMethyl-Histone H3-K27 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; TriMethyl-Histone H3-K27
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human histone H3 (NP_003520.1).
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IP, IF/ICC, IHC-P, ChIP, ChIP-seq, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag