

Trimethyl-Histone H3-K36 (K36Me3) Rabbit Mab, Q16695/P68431 - Wb, Chip-Seq, Ihc-P

Item Number: CM0000295



Introduction

TriMethyl-Histone H3-K36 rabbit mAb (Q16695/P68431) for WB, ChIP-seq, IHC-P. Detects K36me3 modification on histone H3 across human, mouse, rat. Suitable for epigenetics and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	TriMethyl-Histone H3-K36 Rabbit mAb
Remarks/Alias	H3.4; H3/g; H3FT; H3t; HIST3H3; Histone H3; HIST1H3A; TriMethyl-Histone H3-K36
Species	Human
GeneID (human)	8290/ 8350
GeneID	8290/8350
Immunogen	Synthetic peptide A synthetic trimethylated peptide around K36 of human Histone H3 (NP_003520.1).
Source	Rabbit
Category	Methylated Antibodies
Application	WB, DB, IHC-P, IP, ELISA, CHIP, ChIP-seq
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag