

Trimethyl-Histone H3-K36 Rabbit Polyclonal Antibody For Epigenetic Research (Q16695/P68431)

Item Number: CM0000267



Introduction

Rabbit polyclonal antibody targeting trimethylated Lys36 of human histone H3.1t and H3.1 (Q16695/P68431). Validated for WB, ChIP, IF/ICC, and ELISA with broad species reactivity. Ideal for epigenetic and chromatin research.

[Learn More](#)

Parameter	Detail
Product Name	TriMethyl-Histone H3-K36 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; TriMethyl-Histone H3-K36
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human histone H3 (NP_003520.1).
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IF/ICC, IHC-P, DB, ChIP, ChIP-seq, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag