

Trimethyl-Histone H3-K4 (H3K4Me3) Rabbit Polyclonal Antibody - Q16695/P68431 For Chip-Seq & Wb

Item Number: CM0000396



Introduction

Rabbit polyclonal antibody recognizing trimethylated Lys4 of histone H3 (Q16695/P68431). Validated for WB, IHC-P, ChIP-seq, ELISA. Suitable for epigenetic biomarker research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	TriMethyl-Histone H3-K4 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; TriMethyl-Histone H3-K4
Species	Human
GeneID (Human)	H3
GeneID	8290/8350
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human histone H3 (NP_003520.1).
Source	Rabbit
Category	Methylated Antibodies
Application	WB, DB, IHC-P, ELISA, ChIP, ChIP-seq
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag