


Trimethyl-Histone H3-K4 (Q16695 / P68431) Rabbit Monoclonal Antibody For Epigenetic Research & Chip, Wb, Elisa

Item Number: CM0000312



Introduction

Rabbit monoclonal antibody against trimethylated lysine 4 of human Histone H3 (Q16695/P68431). Validated for WB, DB, IHC-P, IF/ICC, ELISA, ChIP. Cross-reacts with Human, Mouse, Rat. Ideal for epigenetic studies and IVD assay development.

[Learn More](#)

Parameter	Description
Product Name	TriMethyl-Histone H3-K4 Rabbit mAb
Remarks/Alias	H3.4; H3/g; H3FT; H3t; HIST3H3; Histone H3; HIST1H3A
Species	Human
GeneID (human)	8290/ 8350
GeneID	8290/8350
Immunogen	Synthetic peptide; A synthetic trimethylated peptide around K4 of human Histone H3 (NP_003520.1)
Source	Rabbit
Category	Methylated Antibodies
Application	WB, DB, IHC-P, IF/ICC, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag