

## Trimethyl-Histone H3-K4 Rabbit Mab - Q16695/P68431 - Wb, IHC-P, If, Chip Antibody

Item Number: CM0000330



### Introduction

Rabbit monoclonal antibody recognizing trimethylated histone H3 at lysine 4 (H3K4me3). Suitable for WB, IHC-P, IF/ICC, ELISA, and ChIP. Predicted cross-reactivity with human, mouse, rat. 15kDa protein weight.

[Learn More](#)

| Parameter        | Value                                                                         |
|------------------|-------------------------------------------------------------------------------|
| Product Name     | TriMethyl-Histone H3-K4 Rabbit mAb                                            |
| Remarks/Alias    | H3.4; H3/g; H3FT; H3t; HIST3H3; Histone H3; HIST1H3A; TriMethyl-Histone H3-K4 |
| Species          | Human                                                                         |
| GeneID (human)   | 8290/ 8350                                                                    |
| GeneID           | 8290/8350                                                                     |
| Immunogen        | A synthetic trimethylated peptide around K4 of human Histone H3 (NP_003520.1) |
| Source           | Rabbit                                                                        |
| Category         | Methylated Antibodies                                                         |
| Application      | WB, IHC-P, IF/ICC, ELISA, ChIP                                                |
| Cross Reactivity | Human, Mouse, Rat, Other (Wide Range Predicted)                               |
| SWISS            | Q16695/P68431                                                                 |
| Protein Weight   | 15kDa                                                                         |
| Shipping         | Ice bag                                                                       |