

Trimethyl-Histone H3-K56 Rabbit Polyclonal Antibody (Q16695/P68431) For Wb, If/Icc, Elisa, Chip-Seq

Item Number: CM0000296



Introduction

High-specificity TriMethyl-Histone H3-K56 rabbit polyclonal antibody recognizing trimethylated lysine 56 of human H3.1/H3.1t. Validated for WB, IF/ICC, ELISA, and ChIP-seq; cross-reacts with human, mouse, rat, and others.

[Learn More](#)

Parameter	Value
Product Name	TriMethyl-Histone H3-K56 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; TriMethyl-Histone H3-K56
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	A synthetic trimethylated peptide around K56 of human Histone H3 (NP_003520.1)
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IF/ICC, ELISA, ChIP-seq
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag