

Trimethyl-Histone H3-K64 (Q16695/P68431) Rabbit Polyclonal Antibody (Wb, Ihc-P, If/Icc, Elisa)

Item Number: CM0000298



Introduction

TriMethyl-Histone H3-K64 rabbit polyclonal antibody for Western blot, IHC-paraffin, IF/ICC, and ELISA. Predicted to recognize trimethylated lysine 64 of histone H3 from human, mouse, rat, and other species. Ideal for epigenetic and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	TriMethyl-Histone H3-K64 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; TriMethyl-Histone H3-K64
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	Synthetic peptide, A synthetic trimethylated peptide around K56 of human TriMethyl-Histone H3-K64 (NP_003520.1).
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag