

Trimethyl-Histone H3-K79 Rabbit Polyclonal Antibody (Q16695 / P68431) For Wb, Chip, Epigenetics

Item Number: CM0000327



Introduction

Rabbit polyclonal antibody specific to trimethylated Histone H3 at lysine 79 (H3K79me3). Suitable for WB, DB, IHC-P, ELISA, ChIP. Cross-reacts with human, mouse, rat. Key tool for chromatin and epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	TriMethyl-Histone H3-K79 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; TriMethyl-Histone H3-K79
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	A synthetic trimethylated peptide around K79 of human histone H3 (NP_003520.1)
Source	Rabbit
Category	Methylated Antibodies
Application	WB, DB, IHC-P, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag