


## Trimethyl-Histone H4-K20 Rabbit Polyclonal Antibody For Wb/Ihc/If (P62805)

Item Number: CM0000257



### Introduction

Explore CamelBio's TriMethyl-Histone H4-K20 Rabbit pAb (P62805) for WB, IHC-P, IF/ICC, DB, ELISA. Specific for human, mouse, rat, and more. Ideal for epigenetic research.

[Learn More](#)

| Parameter        | Value                                                                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | TriMethyl-Histone H4-K20 Rabbit pAb                                                                                                                                |
| Remarks/Alias    | H4; H4/n; H4C1; H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4F2; H4FN; FO108; H4-16; H4C11; H4C12; H4C13; H4C15; H4C16; HIST2H4; HIST2H4A; TriMethyl-Histone H4-K20 |
| Species          | Human                                                                                                                                                              |
| GeneID (human)   | H4                                                                                                                                                                 |
| GeneID           | 8359                                                                                                                                                               |
| Immunogen        | Synthetic peptide; A synthetic trimethylated peptide around K20 of human histone H4 (NP_003529.1).                                                                 |
| Source           | Rabbit                                                                                                                                                             |
| Category         | Methylated Antibodies                                                                                                                                              |
| Application      | WB, DB, IHC-P, IF/ICC, ELISA                                                                                                                                       |
| Cross Reactivity | Human, Mouse, Rat, Other (Wide Range Predicted)                                                                                                                    |
| SWISS            | P62805                                                                                                                                                             |
| Protein Weight   | 11kDa                                                                                                                                                              |
| Shipping         | Ice bag                                                                                                                                                            |