

Triobp Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9H2D6

Item Number: CM0030334

Introduction

TRIOBP polyclonal antibody (Human, Mouse) validated for WB, IHC-P, ELISA. Essential for actin cytoskeleton organization and stereocilia rootlet formation. Ideal for diagnostic and research applications.

[Learn More](#)

Parameter	Value
Product Name	TRIOBP Rabbit pAb
Remarks/Alias	TARA; TAP68; DFNB28; dj37E16.4; HRIHFB2122; TRIOBP
Species	Human
GeneID (Human)	11078
GeneID	11078
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 180-310 of human TRIOBP (NP_008963.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9H2D6
Protein Weight	261kDa
Shipping	Ice bag