


Trka Rabbit Monoclonal Antibody For Wb, Elisa - P04629

Item Number: CM0005670



Introduction

Rabbit monoclonal TrkA antibody for WB and ELISA. Detects human, mouse, and rat TrkA/NTRK1 receptor tyrosine kinase. Ideal for neurotrophic signaling and oncogenic research.

[Learn More](#)

Parameter	Value
Product Name	TrkA Rabbit mAb
Remarks/Alias	MTC; TRK; TRK1; TRKA; Trk-A; p140-TrkA; TrkA
Species	Human
GeneID (Human)	4914
GeneID	4914
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 697-796 of human TrkA (P04629).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P04629
Protein Weight	87kDa
Shipping	Ice bag