


## Trka Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P04629

Item Number: CM0020233



### Introduction

High affinity NGF receptor (TrkA) rabbit polyclonal antibody validated for WB, IHC-P, and ELISA. Reacts with human, mouse, rat. Ideal for neuroscience and oncology research applications. UniProt: P04629.

[Learn More](#)

| Parameter        | Value                                                                                                           |
|------------------|-----------------------------------------------------------------------------------------------------------------|
| Product Name     | TrkA Rabbit pAb                                                                                                 |
| Remarks/Alias    | MTC; TRK; TRK1; TRKA; Trk-A; p140-TrkA; TrkA                                                                    |
| Species          | Human                                                                                                           |
| GeneID (Human)   | 4914                                                                                                            |
| GeneID           | 4914                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human TrkA (NP_002520.2) |
| Source           | Rabbit                                                                                                          |
| Category         | Polyclonal Antibodies                                                                                           |
| Application      | WB, IHC-P, ELISA                                                                                                |
| Cross Reactivity | Human, Mouse, Rat                                                                                               |
| SWISS            | P04629                                                                                                          |
| Protein Weight   | 87kDa                                                                                                           |
| Shipping         | Ice bag                                                                                                         |