


Trka Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P04629

Item Number: CM0013346



Introduction

TrkA Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Targets high-affinity NGF receptor NTRK1. Ideal for neurobiology, oncology research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	TrkA Rabbit pAb
Remarks/Alias	MTC; TRK; TRK1; TRKA; Trk-A; p140-TrkA; TrkA
Target Species	Human
Gene ID	4914
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 700-796 of human TrkA (NP_002520.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04629
Protein Weight	87kDa
Shipping	Ice bag