


## Trkb Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q16620

Item Number: CM0002624



### Introduction

Recombinant rabbit monoclonal antibody against human TrkB (UniProt Q16620). Validated for WB, IHC-P, ELISA in human, mouse, rat. Ideal for neuroscience research on neurotrophin signaling, differentiation, and synaptic plasticity.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | TrkB Rabbit mAb                                                                                                     |
| Remarks/Alias    | OBHD; TRKB; DEE58; trk-B; EIEE58; GP145-TrkB; TrkB                                                                  |
| Species          | Human                                                                                                               |
| GeneID (Human)   | 4915                                                                                                                |
| GeneID           | 4915                                                                                                                |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 32-430 of human TrkB (NP_001018074.1) |
| Source           | Rabbit                                                                                                              |
| Category         | Monoclonal Antibodies                                                                                               |
| Application      | WB, IHC-P, ELISA                                                                                                    |
| Cross Reactivity | Human, Mouse, Rat                                                                                                   |
| SWISS            | Q16620                                                                                                              |
| Protein Weight   | 92kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |