


Trkb Rabbit Pab For Wb, If/Icc, Ihc-P, Elisa - Q16620

Item Number: CM0020768



Introduction

TrkB Rabbit pAb is a polyclonal antibody raised against human TrkB (aa 32-430). Validated for WB, IF/ICC, IHC-P, ELISA. Reacts with human, mouse, rat. Suitable for neurobiology and receptor signaling research.

[Learn More](#)

Parameter	Value
Product Name	TrkB Rabbit pAb
Remarks/Alias	OBHD; TRKB; DEE58; trk-B; EIEE58; GP145-TrkB; TrkB
Species	Human
GeneID (Human)	4915
GeneID	4915
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16620
Protein Weight	92 kDa
Shipping	Ice bag