


Trpv6 Rabbit Pab - Q9H1D0

Item Number: CM0015632



Introduction

This TRPV6 Rabbit pAb detects human, mouse, and rat TRPV6, a calcium-selective channel critical for intestinal Ca^{2+} uptake and homeostasis. Suitable for WB, IHC-P, and ELISA.

[Learn More](#)

Parameter	Value
Product Name	TRPV6 Rabbit pAb
Remarks/Alias	CAT1; CATL; ZFAB; ECAC2; ABP/ZF; HRPTTN; LP6728; HSA277909; TRPV6
Species	Human
GeneID (Human)	55503
GeneID	55503
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human TRPV6 (NP_061116.5)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H1D0
Protein Weight	87kDa
Shipping	Ice bag