

## Trub1 Rabbit Pab - Q8Wwh5

Item Number: CM0022368



### Introduction

TRUB1 Rabbit pAb is a polyclonal antibody for Western Blot (WB) and ELISA, detecting human, mouse, and rat TRUB1, a pseudouridylate synthase involved in mRNA and tRNA pseudouridylation.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | TRUB1 Rabbit pAb                                                                                                   |
| Remarks/Alias    | PUS4; TRUB1                                                                                                        |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 142940                                                                                                             |
| GeneID           | 142940                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 100-280 of human TRUB1 (NP_631908.1) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | Q8WWH5                                                                                                             |
| Protein Weight   | 37kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |