

Tsg101/Vps23 Rabbit Mab For Wb, Elisa - Q99816

Item Number: CM0010778



Introduction

High-specificity rabbit monoclonal antibody against TSG101 (ESCRT-I complex subunit), validated for WB and ELISA. Reacts with human, mouse, and rat samples. Ideal for vesicular trafficking and cancer research.

[Learn More](#)

Parameter	Value
Product Name	TSG101/VPS23 Rabbit mAb
Remarks/Alias	TSG10; VPS23
Species	Human
GeneID (Human)	7251
GeneID	7251
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99816
Protein Weight	44kDa
Shipping	Ice bag