


Tspan15 Rabbit Pab For Wb And Elisa - O95858

Item Number: CM0013763



Introduction

Rabbit polyclonal antibody targeting Tetraspanin-15 (TSPAN15), validated for WB and ELISA. Suitable for human, mouse, and rat samples. Ideal for studying ADAM10 trafficking, Notch signaling, and cell adhesion.

[Learn More](#)

Parameter	Value
Product Name	TSPAN15 Rabbit pAb
Remarks/Alias	NET7; NET-7; TM4SF15; 2700063A19Rik; TSPAN15
Species	Human
GeneID (Human)	23555
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 110-220 of human TSPAN15 (NP_036471.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95858
Protein Weight	33kDa
Shipping	Ice bag