


Tsr1 Rabbit Polyclonal Antibody For Wb - Q2NI82

Item Number: CM0025899



Introduction

Rabbit polyclonal antibody targeting human TSR1, suited for WB and ELISA. Detects endogenous TSR1 in human and mouse samples; essential tool for ribosome biogenesis research. Shipping: Ice bag.

[Learn More](#)

Parameter	Value
Product Name	TSR1 Rabbit pAb
Remarks/Alias	TSR1
Species	Human
GeneID (Human)	55720
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 455-804 of human TSR1 (NP_060598.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q2NL82
Protein Weight	92kDa
Shipping	Ice bag