

Ttk Rabbit Pab - Ttk

Item Number: CM0019887

Introduction

TTK Rabbit pAb is a polyclonal antibody raised against a recombinant fragment of human TTK protein. Suitable for WB and ELISA applications, it exhibits cross-reactivity with human, mouse, and rat samples, with an expected band size of 96-97 kDa.

[Learn More](#)

Parameter	Value
Product Name	TTK Rabbit pAb
Alias	ESK; PYT; CT96; MPH1; MPS1; MPS1L1; TTK
Species	Human
Gene ID	7272
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 536-789 of human TTK (NP_003309.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
Protein Weight	96kDa/97kDa
Shipping	Ice bag