

Tusc2 Rabbit Pab - 075896

Item Number: CM0030646

Introduction

TUSC2 Rabbit polyclonal antibody (pAb) targeting the tumor suppressor candidate 2 (FUS1) protein. Suitable for Western blot, immunohistochemistry, and ELISA applications with cross-reactivity to human, mouse, and rat species. A reliable tool for cancer biology research.

[Learn More](#)

Parameter	Value
Product Name	TUSC2 Rabbit pAb
Remarks/Alias	PAP; FUS1; PDAP2; C3orf11; TUSC2
Species	Human
GeneID (Human)	11334
GeneID	11334
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human TUSC2 (NP_009206.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075896
Protein Weight	12kDa
Shipping	Ice bag