

Twist2 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Wvj9

Item Number: CM0027582

Introduction

Rabbit polyclonal antibody against human TWIST2 (bHLHa39, Dermo-1), validated for WB and ELISA, cross-reacts with mouse and rat. Suitable for studying transcriptional regulation, osteogenesis, and inflammatory responses.

[Learn More](#)

Parameter	Value
Product Name	Twist2 Rabbit pAb
Remarks/Alias	AMS; FFDD3; BBRsAY; Dermo1; SETLSS; bHLHa39; Twist2
Species	Human
Gene ID (Human)	117581
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Twist2 (NP_476527.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8WVJ9
Protein Weight	18kDa
Shipping	Ice bag