

Tyro3 Rabbit Polyclonal Antibody For Wb, Elisa - Q06418

Item Number: CM0017384



Introduction

Tyro3 Rabbit Polyclonal Antibody validated for WB and ELISA. Detects human Tyro3, cross-reacts with mouse and rat. Ideal for studying Tyro3-mediated signaling, neuroprotection, and viral entry mechanisms in research applications.

[Learn More](#)

Parameter	Value
Product Name	Tyro3 Rabbit pAb
Remarks/Alias	BYK; Dtk; RSE; Rek; Sky; Tif; Etk-2; Tyro3
Species	Human
Gene ID (Human)	7301
Gene ID	7301
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 290-430 of human Tyro3 (NP_006284.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q06418
Protein Weight	97kDa
Shipping	Ice bag