

Ube2D1 Rabbit Pab For Wb, Ihc-P, Elisa - P51668

Item Number: CM0024097



Introduction

Rabbit polyclonal antibody against UBE2D1 (UBCH5/UBC4/5) for Western blot, IHC-P, and ELISA. Recognizes human, mouse, and rat proteins. Useful for studying ubiquitin conjugation and protein degradation pathways.

[Learn More](#)

Parameter	Value
Product Name	UBE2D1 Rabbit pAb
Remarks/Alias	SFT; UBCH5; UBC4/5; UBCH5A; E2(17)KB1; UBE2D1
Species	Human
Gene ID (Human)	7321
Gene ID	7321
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-147 of human UBE2D1 (NP_003329.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51668
Protein Weight	17kDa
Shipping	Ice bag