

Ube2F Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q969M7

Item Number: CM0025491



Introduction

Rabbit polyclonal antibody targeting UBE2F (NEDD8-conjugating enzyme 2). Validated for Western blot, IHC-P, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Suitable for studying neddylation and cullin modification pathways.

[Learn More](#)

Parameter	Value
Product Name	UBE2F Rabbit pAb
Remarks/Alias	NCE2; UBE2F
Species	Human
GeneID	140739
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-185 of human UBE2F (NP_542409.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q969M7
Protein Weight	21kDa
Shipping	Ice bag