


Ube2I Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P63279

Item Number: CM0008970



Introduction

Reliable UBE2I (UBC9) rabbit monoclonal antibody for Western blot, IHC-P, and ELISA. Detects human, mouse, and rat SUMO-conjugating enzyme UBC9, a critical regulator of sumoylation, nuclear architecture, and chromosome segregation.

[Learn More](#)

Parameter	Value
Product Name	UBE2I Rabbit mAb
Remarks/Alias	P18; UBC9; C358B7.1; UBE2I
Species	Human
GeneID (Human)	7329
GeneID	7329
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human UBE2I (P63279)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P63279
Protein Weight	18kDa
Shipping	Ice bag