


Ube2l Rabbit Pab For Wb - P63279

Item Number: CM0020692



Introduction

This rabbit polyclonal antibody against UBE2L (UBC9) is validated for WB, IHC-P, IF/ICC, and ELISA in human, mouse, and rat. Explore SUMOylation and chromosome segregation with high specificity.

[Learn More](#)

Parameter	Value
Product Name	UBE2L Rabbit pAb
Remarks/Alias	P18; UBC9; C358B7.1; UBE2L
Species	Human
GeneID (Human)	7329
GeneID	7329
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-158 of human UBE2L (NP_003336.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P63279
Protein Weight	18kDa
Shipping	Ice bag