

Ube2K Rabbit Polyclonal Antibody For Wb, Elisa - P61086

Item Number: CM0013829



Introduction

High-affinity rabbit polyclonal antibody against human UBE2K (E2-25K, HIP2). Validated for Western blot and ELISA, it is an essential tool for ubiquitin-proteasome research, Huntington's disease studies, and ER-associated degradation analysis.

[Learn More](#)

Parameter	Value
Product Name	UBE2K Rabbit pAb
Remarks/Alias	LIG; HIP2; HYPG; UBC1; E2-25K; UBE2K
Species	Human
Gene ID (Human)	3093
Gene ID	3093
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human UBE2K (NP_005330.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P61086
Protein Weight	22kDa
Shipping	Ice bag