


Ube2W Rabbit Polyclonal Antibody For Wb, Elisa - Q96B02

Item Number: CM0022764



Introduction

This UBE2W rabbit polyclonal antibody detects human and mouse UBE2W (UBC-16), validated for Western blot and ELISA. UBE2W is an E2 ubiquitin-conjugating enzyme involved in N-terminal monoubiquitination and DNA repair.

[Learn More](#)

Parameter	Value
Product Name	UBE2W Rabbit pAb
Remarks/Alias	UBC16; UBC-16; UBE2W
Species	Human
GeneID (human)	55284
GeneID	55284
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-110 of human UBE2W (NP_001001481.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96B02
Protein Weight	17kDa
Shipping	Ice bag