


Ube3A Rabbit Pab - Q05086

Item Number: CM0024656



Introduction

UBE3A rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA; detects a 101 kDa band; cross-reacts with human, mouse, and rat. Ideal for studying ubiquitination, circadian rhythms, synaptic plasticity, and viral oncogenesis.

[Learn More](#)

Parameter	Value
Product Name	UBE3A Rabbit pAb
Remarks/Alias	AS; ANCR; PIX1; E6-AP; HPVE6A; EPVE6AP; UBE3A
Species	Human
GeneID (Human)	7337
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-590 of human UBE3A (NP_570853.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q05086
Protein Weight	101kDa
Shipping	Ice bag