


Ube4A Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q14139

Item Number: CM0028377



Introduction

Rabbit polyclonal antibody against human UBE4A (ubiquitin conjugation factor E4 A) for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for ubiquitin pathway research.

[Learn More](#)

Parameter	Value
Product Name	UBE4A Rabbit pAb
Remarks/Alias	E4; UFD2; UBOX2; NEDHMS; UBE4A
Species	Human
Gene ID (Human)	9354
Gene ID	9354
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 797-1066 of human UBE4A (NP_001191006.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14139
Protein Weight	123kDa
Shipping	Ice bag