

## Ubiquitin Rabbit Pab - P0Cg47

Item Number: CM0015192



### Introduction

Ubiquitin rabbit polyclonal antibody for WB, IHC-P, ELISA. Reacts with human, mouse, and rat ubiquitin. Ideal for studying proteasomal degradation, DNA repair, and signal transduction.

[Learn More](#)

| Parameter        | Value                                                                                                                |
|------------------|----------------------------------------------------------------------------------------------------------------------|
| Product Name     | Ubiquitin Rabbit pAb                                                                                                 |
| Remarks/Alias    | HEL-S-50; Ubiquitin                                                                                                  |
| Species          | Human                                                                                                                |
| GeneID (human)   | 7314                                                                                                                 |
| GeneID           | 7314                                                                                                                 |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-228 of human Ubiquitin (NP_061828.1) |
| Source           | Rabbit                                                                                                               |
| Category         | Polyclonal Antibodies                                                                                                |
| Application      | WB, IHC-P, ELISA                                                                                                     |
| Cross Reactivity | Human, Mouse, Rat                                                                                                    |
| SWISS            | P0CG47                                                                                                               |
| Protein Weight   | 26kDa                                                                                                                |
| Shipping         | Ice bag                                                                                                              |