


Ubiquitin Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P0CG47

Item Number: CM0022900



Introduction

High-quality Ubiquitin Rabbit Polyclonal Antibody suitable for WB, IHC-P, IF/ICC, and ELISA applications. Detects ubiquitin in human, mouse, rat, and other species. Ideal for protein degradation and cell signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Ubiquitin Rabbit pAb
Remarks/Alias	HEL-S-50; Ubiquitin
GeneID (Human)	7314
GeneID	7314
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-228 of human Ubiquitin (NP_061828.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P0CG47
Protein Weight	26kDa
Shipping	Ice bag