


Ubiquityl-Histone H2A-K119 (P04908/P0C0S8) Rabbit Polyclonal Antibody For Wb And Elisa

Item Number: CM0000075



Introduction

Rabbit polyclonal antibody targeting ubiquitylated histone H2A at lysine 119 (UniProt P04908/P0C0S8). Suitable for Western blot and ELISA, with wide species reactivity (Human, Mouse, Rat, and predicted). Ideal for epigenetics and chromatin research in IVD development.

[Learn More](#)

Parameter	Value
Product Name	Ubiquityl-Histone H2A-K119 Rabbit pAb
Remarks/Alias	HIST1H2AB; Ubiquityl-Histone H2A-K119
Species	Human
GeneID (human)	H2A
GeneID	3012/8329
Immunogen	Synthetic Peptide: A synthetic ubiquitinated peptide around K119 of human HIST1H2AB (NP_734466.1).
Source	Rabbit
Category	Other Modified Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P04908/P0C0S8
Protein Weight	14kDa
Shipping	Ice bag