


Uch37/Uchl5 Rabbit Mab For Wb And Elisa - Q9Y5K5

Item Number: CM0007414



Introduction

Rabbit monoclonal antibody against human UCH37/UCHL5, validated for WB and ELISA. Recognizes human, mouse, and rat UCHL5, a proteasome-associated deubiquitinase involved in protein degradation and DNA repair.

[Learn More](#)

Parameter	Value
Product Name	UCH37/UCHL5 Rabbit mAb
Remarks/Alias	UCHL5; CGI-70; INO80R; UCH-L5; UCH37; ubiquitin C-terminal hydrolase L5; UCH37/UCHL5
Species	Human
Gene ID	51377
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-329 of human UCH37/UCHL5 (NP_057068.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5K5
Protein Weight	36kDa/37kDa
Shipping	Ice bag